

New! Low Temperature Coolers for H/D Exchange

- Temperature range between ambient and 2 °C.
- Temperature stability to 0.2 °C
- 0 °C possible
- Columns enter into and exit from the cooler through fittings
- Advanced thermoelectric cooling technology

Model number:


PST-CC column cooler for multiple columns. Cooling volume: 10.2 cm x 8 cm x 1.2 cm

PST-CC-01 Butterfly (flexible) column cooler for transfer tubing and columns, 13 cm x 10.5 cm

Note: Outside dimensions are larger. Dimensions of heat sink and fan are not included.


PST-CHC-LT Column cooler controller with power supply

Please call to inquire about custom installation on nanospray sources



Also available:

PST-CH-LT : column heater used inside a cold box



- 5 cm x 5 cm x ~1 cm outside dimensions
- Accommodates columns of various lengths by coiling the column
- For sample prep capillary as well as separation column

© 2012 Phoenix S&T | All rights reserved